

PETRISWISS MINI PS20

The new PETRISWISS MINI PS20 is a fully automatic and very precise Petri dish filling system, "made in Switzerland", of the most recent generation, based on state-of-the-art findings, featuring a broad range of possible applications and with a sensational price-performance ratio.

FIELDS OF APPLICATION

Laboratories, research and development, the chemical industry, the pharmaceutical industry, hospitals, the food industry and universities etc.

POSSIBLE APPLICATIONS

DISH FEED	Feed turret with 20 Petri dishes
SINGLE DOSING	Programming medium dosing quantity per Petri dish
STACKING	Stacking quantity (1-20) of the filled Petri dishes
DISH DISCHARGE	Automatic discharge

INCLUDING USEFUL FEATURES !

All program data, parameters and operating states are displayed clearly and understandably on a graphic display. The unit is very operator-friendly. The unit is operated using a state-of-the-art touch screen.

- Graphic display with backlighting for perfect readoff.
- Easily understandable, user-friendly operation with touch screen.
- Operator prompting in German or English, freely selectable by the user.
- Fully automatic filling and stacking of the Petri dishes.
- High, trouble-free Petri dish throughput.
- Dosing quantities, speed, pause time and AntiDrop function programmable.
- Stack function.
- Automatic self-test initialisation and run.
- Self-monitoring of all device functions (motors and switches).
- Easy care and servicing of the individual device components.

