



LF.8.201.000 standard

The LabFlame basic can be operated using the included foot pedal (with cable) or the button function.

Various functions are available: flexible start-stop function, foot pedal control, and button control.



LF. 8.200.000

The LabFlame IR is the model that offers the most customization options and is our most versatile model with additional safety features:

IR Features:

- Display with self-explanatory icons
- Acoustic warning signals
- IR sensor for activation through hand movements
- Gas consumption indicator

Safety Features:

- No accidental activation with double-click
 IR sensor
- Zero-pressure shutdown
- Cooling reminder

Customizations:

- Up to two users
- Individual setting of double-click function
- Individual temperature adjustment
- Individual cooling time
- Individual burning time



LF. 8.203.000

Wireless and secure. With the radio frequency (RF) foot pedal, the LabFlame RF can be operated without cables.

Other than that, the LabFlame RF has the same functions as the standard model.



LF. 8.204.000

Equipped with two (AA) batteries, the LabFlame + transforms into the most mobile device, ready for flexible use.

With a single battery charge, continuous operation of up to 9 hours is possible (2000 inoculation loops).

Other than that, the LabFlame + has the same functions as the IR model.



BioTool AG Industrie Neuhof 9 3422 Kirchberg Switzerland

+41 34 423 50 60 info@biotoolswiss.com www.biotoolswiss.com